

https://community.nanoporetech.com/protocols/genomic-dna-sqk-nsk007/v/gde_9002_v7_rev16may2016-449/dna-repair

- **Perform FFPE DNA repair treatment using NEB M6630.**

Reagent	Volume
1-1.5 μ g fragmented** DNA	45 μ l
Nuclease-free water	8.5 μ l
FFPE Repair Buffer	6.5 μ l
FFPE Repair Mix	2 μ l
Total	62 μ l

- **Mix gently by inversion and spin down.**
- **Incubate the reaction for 15 minutes at 20 °C.**
- **Resuspend AMPure XP beads by vortexing.**
- **Add 62 μ l of the resuspended beads to the End-prep reaction and mix gently by pipetting.**
- **Incubate on a rotator mixer (Hula mixer) for 5 minutes at *room temperature*.**
- **Rotator mixer**

The agitation is gentle and the contents of the tube may not appear to move at all. This allows the DNA to bind to the beads.

- **Prepare 500 μ l of fresh 70% ethanol in nuclease-free water.**
- **Spin down the sample and *pellet on a magnet*. Leaving the tube on the magnet, pipette off the supernatant.**
- **Keep on magnet, wash beads with 200 μ l of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard. Repeat.**
- **Spin down and replace on magnet to collect and pipette off any residual 70% ethanol. Briefly allow to dry.**
- **Remove the tube from the magnetic rack and resuspend pellet in 46 μ l nuclease-free water. Incubate for 2 minutes at *room temperature*.**
- ***Pellet beads* on magnet until the eluate is clear and colourless.**
- **Remove and retain 46 μ l of eluate in a clean 1.5 ml Eppendorf DNA LoBind tube.**