

Lambda control experiment for the MinION™ device (1/2)

Flow Cell Number

DNA Samples

Before start checklist

- ☐ Vortexer and Hula mixer
- ☐ Pipettes P1000, P200, P100, P20, P10 and P2
- ☐ Tips ready: 40 x 1000 µl, 30 x 200 µl, 20 x 10 µl
- ☐ Covaris g-TUBE
- ☐ AMPure beads resuspended and at RT
- ☐ MyOne C1 Streptavidin beads resuspended and at RT

- ☐ Magnet for bead separation
- ☐ Approx 20 DNA LoBind Eppendorf tubes
- ☐ LMD / BBB / ELB / DCS / HPT / RBF1 on ice
- ☐ HPA / AMX in freezer until needed
- ☐ Freshly prepared 70% EtOH
- ☐ Nuclease free water (NFW)
- ☐ Notebook sleeptimer / update off
- ☐ 10mM Tris-HCl pH 8.5
- ☐ Microfuge



- ☐ Nanopore sequencing kit SQK-NSK007
- ☐ Update checks for MinKNOW™ and Metrichor™ completed
- ☐ Run platform QC prior to library prep
- ☐ NEB Blunt / TA Ligase Master Mix
- ☐ 0.2 ml thin-walled PCR tubes
- ☐ Heat block at 37°C capable of taking 1.5ml tubes
- ☐ Thermal cycler AT 20°C & 65°C
- ☐ NEBNext Ultra II End-repair / dA-tailing Module

MASSFLOW	INSTRUCTIONS	NOTES / OBSERVATIONS	TIME / DATE
	Lambda control experiment 1 µg Genomic DNA in 46 µl (20 µl LMD and 26 µl NFW) Shear in a Covaris g-TUBE Eppendorf 5424; 6000 rpm for 8 kb fragments. 2 x 1 minute (invert tube to collect) Keep processing time under 15 minutes	A PreCR step is not required when using the Lambda DNA contained in the kit as it is high quality DNA, however it is recommended with other gDNA samples	
	1 µl Agilent Bioanalyzer / Gel analysis		
	Add from NEBNext Ultra II End-Repair / dA-tailing Module <input type="checkbox"/> 7 µl Ultra II End-Prep buffer <input type="checkbox"/> 3 µl Ultra II End-Prep enzyme mix <input type="checkbox"/> 5 µl DNA CS 3.6 kb (DNA CS – positive control) Mix by inversion + spin down Incubate for at 20 °C for 5 minutes and 65 °C for 5 minutes		
	<input type="checkbox"/> Add 60 µl resuspended AMPure XP beads at RT <input type="checkbox"/> Incubate on rotator for 5 minutes, spin down and pellet on magnet. Discard the supernatant. <input type="checkbox"/> Keep on magnet, wash 2x with 200 µl fresh 70% EtOH, do not disturb pellet <input type="checkbox"/> Briefly spin down, replace on magnet, pipette off residual wash. Briefly allow to dry <input type="checkbox"/> Resuspend pellet in 31 µl Nuclease-free water, incubate at RT for 2 minutes <input type="checkbox"/> Pellet beads on a magnet, remove eluate of End-Prepped DNA and transfer to fresh DNA LoBind		
	1 µl QuBit fluorimeter – recovery aim about 700 ng of material		
	NEB Blunt / TA Ligase Master Mix. When preparing the ligation reaction, mix by inversion between each addition. <input type="checkbox"/> Check that the Master Mix is clear of any precipitate <input type="checkbox"/> 8 µl NFW <input type="checkbox"/> 10 µl Adapter Mix <input type="checkbox"/> 2 µl HPA <input type="checkbox"/> 50 µl NEB Blunt / TA Master Mix Incubate at RT for 10 minutes		
	<input type="checkbox"/> 1 µl HPT Mix by inversion + spin down Incubate at RT for 10 minutes While incubating prepare MyOne C1 beads for purification of the adapted and tethered DNA		
	Prepare the MyOne C1 beads in DNA LoBind tube <input type="checkbox"/> Vortex MyOne C1 beads to resuspend <input type="checkbox"/> Transfer 50 µl of beads to DNA LoBind tube <input type="checkbox"/> Pellet beads on magnet and remove supernatant <input type="checkbox"/> Wash beads 2x with 100 µl BBB vortexing to resuspend. Pellet on magnet, pipette off and discard supernatant <input type="checkbox"/> Resuspend washed beads in 100 µl BBB to be used for adapter ligation cleanup		
	Purifying the adapted DNA <input type="checkbox"/> Add 100 µl washed MyOne C1 beads to the adapted, tethered DNA reaction on a rotator <input type="checkbox"/> Incubate at RT for 5 minutes		
	<input type="checkbox"/> Place on magnet, pellet and pipette off the supernatant <input type="checkbox"/> Wash beads with library bound 2x with 150 µl BBB, pipette to resuspend. Pellet on magnet pipette off and discard supernatant <input type="checkbox"/> Close the lid, spin down, replace on magnet for 1-2 minutes, open lid and pipette off any residual BBB		
	Elute the adapted library (Pre-sequencing Mix) <input type="checkbox"/> Resuspend pelleted beads in 25 µl ELB and incubate at 37°C for 10 minutes <input type="checkbox"/> Pellet beads on the magnet, remove the eluate and transfer to new DNA LoBind tube		

Burn-In for the MinION™ device (2/2)

Flow Cell Number

DNA Samples

MASSFLOW	INSTRUCTIONS	NOTES / OBSERVATIONS	TIME / DATE
	1 µl QuBit fluorimeter – recovery aim ~250ng ng of material		
Before start checklist <div> <input type="checkbox"/> MinION™ connected to computer with MinION™ Flow Cell <input type="checkbox"/> Run platform QC in parallel to library prep <input type="checkbox"/> MinION with MinION Flow Cell </div> <div> <input type="checkbox"/> Computer setup to run MinKNOW <input type="checkbox"/> Metrichor Agent setup <input type="checkbox"/> Run Name set <input type="checkbox"/> Screws and Heat pad for original MinION </div> <div> <input type="checkbox"/> Pre-sequencing Mix (PSM) at > 4 ng/µl <input type="checkbox"/> PSM and RBF1 on ice <input type="checkbox"/> NFW at RT <input type="checkbox"/> Platform QC completed </div>			
	Prepare the MinION for sequencing protocol This step can be run in parallel with the preparation of the library from genomic DNA to Presequencing Mix <input type="checkbox"/> Assemble the MinION and MinION Flow Cell <input type="checkbox"/> Setup MinKNOW to run the Platform QC – name the run and start the protocol script – NC_Platform_QC_FLO_MIN104.py <input type="checkbox"/> Allow the script to run to completion and the number of active pores are reported		
	Prime the flow cell ready for the library to be loaded Prepare priming buffer <input type="checkbox"/> 500 µl 2x RBF1 <input type="checkbox"/> 500 µl NFW Mix by vortexing and spin down	RBF1	
	Prime the flow cell <input type="checkbox"/> Open the sample port and using a 1000 µl pipette draw back the buffer in the flow cell to make sure there is continuous buffer flow from the sample port across the sensor array. <input type="checkbox"/> Load 500 µl of the priming buffer using a vertical 1000 µl pipette and tip in a continuous flow avoiding introducing air bubbles or disturbing the sensor array. Wait 10 minutes, then repeat the loading and wait another 10 minutes		
	Loading the prepared library <input type="checkbox"/> Mix the Pre-sequencing Mix by inversion 10x, briefly spin down the library <input type="checkbox"/> Add the reagents in the following order <input type="checkbox"/> 75 µl RBF1 kept at RT <input type="checkbox"/> 63 µl NFW kept at RT <input type="checkbox"/> 12 µl Pre-sequencing Mix kept on ice <input type="checkbox"/> Load 150 µl of the library loading mix via the sample port, close the cover and replace the MinION lid.	RBF1	
	Starting the sequencing script in MinKNOW and the workflow in the Metrichor Agent <input type="checkbox"/> Return to MinKNOW, name the run, select the NC_6Hr_Lambda_Burn_In_Run_FLO_MIN104.py and start using the start in the MinKNOW dialogue box <input type="checkbox"/> Open the Metrichor Agent, select the latest version of the Lambda Control Experiment workflow, run the workflow and monitor the workflow using <input type="checkbox"/> the visualisation options in details MinKNOW will report the number of pores available for sequencing before data collection begins. These may differ from those reported in the Platform QC. <input type="checkbox"/> Allow the protocol to proceed until MinKNOW reports Finished Successfully System Ready . Use the Stop in the Control Panel <input type="checkbox"/> to finish the protocol. Quit the Metrichor Agent, close down MinKNOW and disconnect the MinION.		
After sequencing checklist <div> <input type="checkbox"/> Store washed flow cell at 4°C or complete the returns form in the Nanopore Community <input type="checkbox"/> Store MinION at RT <input type="checkbox"/> Return reagents to the freezer </div> <div> <input type="checkbox"/> Navigate to www.metrichor.com to review the full sequencing report </div>			