

Equipment and consumables for Rapid Sequencing Kit



Oxford Nanopore Technologies
Oxford Science Park, Oxford OX4 4GA, UK
support@nanoporetech.com
www.nanoporetech.com

Equipment and consumables for Rapid Sequencing Kit

Lambda control experiment

Equipment	Hula mixer (gentle rotator mixer)
	Microfuge
	Vortex mixer
	Thermal cycler at 30 °C and 75 °C
	P1000 / (P200 / P100) / (P20 / P10) / P2 pipette
	Ice bucket with ice
	Pre-chilled freezer block at -20 °C for 200 µl tubes
	Eppendorf tube marker pen
	Timer
Third party consumables	NEB Blunt/TA Ligase Master Mix (M0367)
Consumables	Nanopore Sequencing Kit
	0.2 ml thin-walled PCR tubes
	Nuclease-free water
	P1000 / (P200 / P100) / (P20 / P10) /P2 pipette tips